



Models	861388 (text display) 861389 (ECG and text display) Supplied with AED, primary battery (1), SMART Pads III (1 set), printed instructions (Setup Guide) and CD-ROM (Administrative Reference)
Waveform	SMART biphasic Truncated Exponential waveform parameters adjust as a function of patient impedance. Adult nominal peak current 32A (150J into a 50 ohm load); pediatric nominal peak current 19A (50J into a 50 ohm load) using optional Infant/Child Key
Shock delivery	Via defibrillator pads placed in the anterior-anterior (Lead II) position for adults; anterior-posterior position for infants and children under 55 lbs (25 kg) or 8 years old
Controls	On/Off button, shock button, option buttons. Auto-On feature, when used with the optional FR3 carry case, enables FR3 to power up when case lid is opened
Indicators	High-resolution color LCD, beeper, voice prompts, tones and chirps, audio speaker, connector socket, ready light, shock button
Advanced Mode	Configurable using optional heartStart Configure software

ECG Display

Screen	LCD color display, 320 x 240 pixels. 7.2 cm x 5.4 cm
Bandwidth	1 Hz to 30 Hz (-3dB), nominal (non-diagnostic)
Monitored Lead	Lead II using anterior-anterior adult pads placement

Physical

Size	6.9 cm x 13.5 cm x 22.1 cm
Weight	1.6 kg with FR3 standard battery installed

Environmental/physical requirements