

ECG RECORDER COMPATIBILITY LIST

Decapus® III

Please contact Quickels Systems AB if you can not find your ECG printer of choice within the list below.



RECORDER CABLE
(QN NO)

CONNECTOR

COMPATABLE ECG RECORDERS

QN D13SC (1.3M)
QN D27SC (2.7M)
QN D40SC (4.0M)

15 D-SUB w. 2 screws

QN D27SC_R (2.7M)
Note: Reversed order
for ECG leads



SCHILLER:

- Cardiovit AT-1
- Cardiovit AT-2
- Cardiovit AT-2 plus
- Cardiovit AT-10 plus
- Cardiovit AT-101
- Cardiovit AT-102
- Cardiovit AT-102 plus
- Cardiovit AT-104 PC
- Cardiovit CS-200
- Cardiovit CS-200 Excellence
- Cardiovit CS-200 Office
- Cardiovit CS-200 Touch

MORTARA:

- ELI 10
- ELI 150c | Adapter QN D752 is needed
- ELI 230 | Adapter QN D752 is needed
- ELI 250c | Adapter QN D752 is needed
- ELI 280 | Adapter QN D752 is needed
- ELI 350 (Stress) | Adapter QN D752 is needed
- ELI 380 | Adapter QN D752 is needed

BTL:

- BTL-08 LC Plus ECG
- BTL-08 LC ECG
- BTL-08 LT Plus ECG
- BTL-08 LT ECG
- BTL-08 MT Plus
- BTL CardioPoint-ECG C600

QUICKELS SYSTEMS AB
Björnnäsvägen 21
SE-114 19 STOCKHOLM
Phone: +46 (0)8 709 49 00
Fax: +46 (0)8 709 00 12
E-mail: info@quickels.com
www.quickels.com

Changes may occur without prior notice.

ver.15.22.4

RECORDER CABLE
(QN NO)

CONNECTOR

COMPATABLE ECG RECORDERS

- BTL CardioPoint-ECG C300
- BTL-08 SD6
- BTL-08 SD3
- BTL-08 SD1

EDAN:

- EDAN DP12 sampling box

PHILIPS:

- Page Writer TC70 | Adapter QN D305-4 is needed
- Page Writer TC50 | Adapter QN D305-4 is needed
- Page Writer TC30 | Adapter QN D305-4 is needed
- Page Writer TC10 | Adapter QN D305-4 is needed

RECORDER CABLE (QN NO)	CONNECTOR	COMPATABLE ECG RECORDERS
---------------------------	-----------	--------------------------

QN D13SI (1.3M)
QN D27SI (2.7M)
QN D40SI (4.0M)

15 D-SUB Lateral grip

SIEMENS:
• Megacart

QN D27SI_R (2.7M)
Note: Reversed
order for ECG leads

CARDIOLEX:
• Lexor

WELCH ALLYN
• Cardio Perfect MD



RECORDER CABLE
(QN NO)

CONNECTOR

COMPATABLE ECG RECORDERS

QN D27WA (2.7M)

REDEL




WELCH ALLYN:

- Cardio Perfect PRO
-



ECG ADAPTERS

Decapus® III

ECG ADAPTER (QN NO)	NAME	DESCRIPTION
QN D305-3	GE adapter series	ECG adapter suitable to connect to ECG printers of type GE MAC 5000.
		
QN D305-4	TC adapter series	ECG adapter suitable to connect to ECG printers of type Philips Page Writer TC Series.
		
QN D752	Mortara adapter	ECG adapter suitable to connect to ECG printers of type Mortara ELI Series.
		

In order to connect Decapus® III to an ECG printer via above mentioned ECG adapters the recorder cable; QN D13SC, QN D27SC or QN D40SC need to be used.

QUICKELS ORDERING SYSTEM (ONT):

Select the recorder cable of choice; QN D13SC, QN D27SC or QN D40SC as well as the ECG adapter when ordering a complete system.